



Product Data Sheet

2304

Rev 2
 Effective date: 20 Dec 2010
 Page 1 of 1
 Authorisation: CO 2021

Product Description

The Browne Washer Disinfectant soil test has been carefully designed to simulate the soiling that occurs during theatre use and performs equivalent to the Edinburgh Test Soil. The Browne Soil Test is supplied in individual pots; water is added to the pot then shaken vigorously to form a paste which is then applied to a test load with the brush provided. After a drying period the test load should be cleaned using the normal procedures and inspected for residual soil. The bright red colour of the Browne Test Soil allows easy identification of areas that have not been properly cleaned.



Specification

The Browne Test Soil has the equivalent performance to the Edinburgh Test Soil described in BS 2745, EN ISO 15883-5 & HTM 2030

Applicable standards

BS2745:1993 specification for Washer-disinfectors for medical purposes.
 EN ISO 15883-5
 HTM2030 Washer-disinfectors.

Product Size

N/A

Packaging

Primary case	10 Bottles & IFU	Size:	245mm x 100mm x 70mm	Weight:	420g Gross
contents:					

Recommended Storage

0 - 5°C 30 – 60 % RH

Shelf Life

12 months from date of manufacture (Store in refrigerator)

Active Components

The dried test soil formula contains protein, lipid and polysaccharide. None of the substances used in the formulation are known carcinogens, nor do they contain any heavy metals. If contact is made wash the skin with soap and water.

If ingested wash out the mouth thoroughly and give plenty of water to drink.

Safe Usage

When reconstituted the Washer Disinfectant Soil contains inert ingredients, colouring dyes, egg proteinaceous materials. As such it does not present a risk to health when used as directed.

Avoid contact with skin and eyes. Do not ingest. If material inadvertently gets into eyes flush with water and seek medical advice.

The product and packaging do not contain latex.

Contains no blood products.

Disposal

Dispose of used/expired soil, bottles and boxes using normal waste disposal.

