

What components of Dietary Fiber does each method measure accurately?


AOAC Method	Megazyme Product range	Target Analytes				Issues	
		Soluble DF		Insoluble DF			
		HMWDF	SDFS (NDO)	SDFP	IDF	Underestimated	Overestimated
985.29	K-TDFR	•				RS ₂ , RS ₃ NDO not measured	RS ₄
991.43				•	•		
2009.01	K-INTDF	•	•			RS ₂ , RS ₄ , FOS	Resistant maltodextrins artifacts
2011.25			•	•	•		
2017.16	K-RINTDF	•	•			NONE	NONE
2022.01			•	•	•		

DF: Dietary Fiber. HMWDF: High Molecular Weight Dietary Fiber. RS: Resistant Starch. IDF: Water insoluble Dietary Fiber. SDFP: Water soluble Dietary Fiber which precipitates in 78% ethanol. SDFS: Water soluble Dietary Fiber that remains soluble in 78% ethanol (= NDO).

• Target analytes

Which AOAC method is most appropriate to measure Dietary Fiber?

What is the Fiber Profile?



DF: Dietary Fiber

RS: Resistant Starch

NDO: Non-digestible Oligosaccharides (=SDFS)

Non-digestible Oligosaccharides	Total Dietary Fiber (TDF)	2017.16	2009.01	TDF	2017.16	2001.03	985.29
Inulin/Fructooligosaccharides (FOS)	Insoluble & Soluble DF	2022.01	2011.25	Insoluble & Soluble DF	2022.01	2001.03*	991.43
Galactooligosaccharides (GOS)	Megazyme Product range	K-RINTDF [†]	K-INTDF	Megazyme Product range	K-RINTDF	Individual Enzymes [#]	K-TDFR
Polydextrose (PDX)							
Resistant Maltodextrins (RMD)							

[†] K-RINTDF provides a more accurate measurement for Resistant Starch

[#] E-AMGDFPD, E-BLAAM & E-BSPRPD

* Modification as per AOAC 991.43 to allow for soluble/insoluble DF determination