

3M Food Safety

3M™ Sample Collection Products



Experience confidence

The complete 3M sample handling system

Collect.

Better testing begins with better samples — and better samples are obtained when you use the best collection tools. 3M™ Sample Collection solutions provide quality products designed specifically for microbial surface sampling. Make environmental testing and product sampling easy and convenient with the full line of 3M Sample Collection products.

- Easy to use for product sampling and environmental testing
- Exceptional convenience and quality
- Increased productivity and efficiency
- Reliable and accurate — time after time
- Improved consistency from technician to technician and plant to plant



Collect.

Sponges and swabs designed specifically for microbial surface sampling.

3M™ Sponge-Sticks*

- Allows sampling without the need to directly handle the sponge, increasing operator efficiency with less risk of contamination
- Easier to reach ceilings, into drains and pipes, and around equipment
- Sturdy sponge snaps free from plastic handle after sampling
- Biocide-free cellulose sponge maintains *Listeria* organism viability within 72 hours of refrigeration**
- Available in dry or pre-hydrated formats, with or without irradiated gloves



3M™ Hydrated Sponges

- Ready-to-use, pre-hydrated sponges for easy and convenient sample collection
- Available in a variety of pre-hydrated dilutions
- Biocide-free cellulose sponge maintains *Listeria* organism viability within 72 hours of refrigeration**
- Durable sponge withstands scrubbing
- Sponge bag with write-on label
- Available with or without irradiated gloves




3M™ Dry-Sponges

- Dry sponge in a sample bag
- Ideal for those who prefer to use their own sampling broths or a broth not available with the 3M™ Hydrated Sponges
- Durable sponge made from biocide-free cellulose to help maintain organism viability



3M™ Extendable Poles

- Designed for use with the 3M™ Sponge Stick for easier sampling of remote and hard-to-reach areas
- Lightweight, sturdy aluminum construction
- Minimizes bending and stretching of plant personnel
- Safely sample ceilings and tanks



*Formerly known as Spongesicle.

**Technical data summary can be provided upon request.

Catalog No.	Description	Qty.
3M™ Sponge-Sticks		
SSL100	3M™ Sponge-Stick, 1.5" x 3" Sponge, dry, 30 oz. bag	100
SSL10BPW	3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Buffered Peptone Water Broth, 30 oz. bag, 5 bags of 20	100
SSL10LET	3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Lethen Broth, 5 bags of 20	100
SSL10LET2G	3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Lethen Broth, 30 oz. bag, 2 gloves, 5 bags of 20	100
SSL10DE	3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL D/E Neutralizing Broth, 30 oz. bag, 5 bags of 20	100
SSL10NB	3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 30 oz. bag, 5 bags of 20	100
SSL10NB2G	3M™ Sponge-Stick, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 30 oz. bag, 2 gloves, 5 bags of 20	100

3M™ Hydrated Sponges		
HS10BPW	3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Buffered Peptone Water Broth, 5 bags of 20	100
HS10BPW2G	3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Buffered Peptone Water Broth, 2 gloves, 5 bags of 20	100
HS10DE2G	3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL D/E Neutralizing Buffer, 2 gloves	100
HS10LET	3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Lethen Broth, 5 bags of 20	100
HS10LET2G	3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Lethen Broth, 2 gloves, 5 bags of 20	100
HS119510NB	3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 50 oz. bag	100
HS10NB	3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 5 bags of 20	100
HS10NB2G	3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 2 gloves, 5 bags of 20	100
HS2410NB2G	3M™ Hydrated Sponge, 1.5" x 3" Sponge w/10mL Neutralizing Buffer, 24 oz. bag, 2 gloves	100

3M™ Dry-Sponges		
BP135DS	3M™ Dry-Sponge w/4" String, 1.5" x 3" Sponge, 18 oz. bag	100
BP133ES	3M™ Dry-Sponge, 18 oz. printed bag	100
BP237SPG	3M™ Dry-Sponge, Sponge w/glove, 18 oz. bag	100

3M™ Extendable Poles		
XSSL10NB	3M™ Extendable Pole Kit, (1) XSL-008 and (1) SSL-10NB	1
XSL008	3M™ Extendable Pole w/Sponge Holder, 8 ft.	1



3M™ Quick Swabs

- Ready-to-use environmental swab quickly delivers 1mL of sample onto a 3M™ Petrifilm™ Plate
- Increased accuracy
- Improved consistency from technician to technician and plant to plant
- Labor savings in dilution mixing and pipetting
- Better productivity and efficiency
- Can be used wet or dry



3M™ Swab-Samplers

- Ideal for sampling hard-to-reach areas with long, pre-hydrated 9cm swab
- Premeasured fill volumes allow for qualitative or quantitative use
- Blue stick is easy to detect if dropped onto production area or floor
- Can stand alone or fit into a standard test tube rack
- Ample headspace for effective vortexing

3M™ Enviro Swabs

- Designed to minimize the risk of sample contamination and improve the detection of bacteria on all types of surfaces
- Large swab head with shaped tip to sample crevices for efficient collection
- Flexible plastic stem allows swabbing of difficult-to-access areas
- Free of substances that inhibit bacteria
- Specially formatted transport medium
- Easily re-sealed for transport



Catalog No.	Description	Qty.
3M™ Quick Swabs		
6432	3M™ Quick Swab	50
6433	3M™ Quick Swab	250
3M™ Swab-Samplers		
RS96010BPW	3M™ Swab-Sampler w/10mL Buffered Peptone Water Broth	100
RS96010DE	3M™ Swab-Sampler w/10mL D/E Neutralizing Broth	100
RS9601LET	3M™ Swab-Sampler w/1mL Lethen Broth	100
RS9604LET	3M™ Swab-Sampler w/4mL Lethen Broth	100
RS9605LET	3M™ Swab-Sampler w/5mL Lethen Broth	100
RS96010LET	3M™ Swab-Sampler w/10mL Lethen Broth	100
RS9604NB	3M™ Swab-Sampler w/4mL Neutralizing Buffer	100
RS96010NB	3M™ Swab-Sampler w/10mL Neutralizing Buffer	100
3M™ Enviro Swabs		
ENVS WB25	3M™ Enviro Swab	25
ENVS WB25D	3M™ Enviro Swab, Dry	25



Global experience. Local support.

3M Food Safety products are at work in more than 60 countries around the globe, to help keep businesses like yours moving forward. Our passion is to work with you to discover new food safety solutions that protect your brand, mitigate your risk and improve your operations. And, every day, our people are there to support you locally, by providing technical help and exceptional service.

For more information, contact your local 3M Food Safety representative or visit 3M.com/foodsafety.



3M Food Safety
3M Center
Building 275-5W-05
St. Paul, MN 55144-1000
USA
1-800-328-6553
www.3M.com/foodsafety

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800-364-3577

Please recycle. Printed in U.S.A. © 3M 2013.
All rights reserved. 70-2011-5044-1

3M and Petrifilm are trademarks of 3M.
Used under license in Canada.